

Product Name

Tracheal Tube (Bougie) - Vented

General Description

The P3 Intubating Introducer consists of a catheter and a single 15mm connector. The catheter has a through lumen design with two distal side ports to allow adequate airflow. The catheter has centimetre markings, which start at 10cm from the distal end, to accurately judge introducer depth within the trachea. The blunt catheter tip is a traumatic.

The 14ch Adult P3 Intubating Introducer has a 4.67 mm outer diameter. The catheter is 750mm long, has a closed end with two side ports. The distal 25 mm portion is angled.

Intended Use

The 14ch Adult P3 Intubating Introducer has been designed for placement of a single lumen tracheal tube whose inner diameter is 6mm or larger. The 14ch Adult P3 Intubating Insufflation Airway Introducer is not intended for use with double lumen tracheal tubes. The removal of the supplied 15mm adapter allows the use of low-pressure oxygen sources, if necessary, during the procedure. The product is intended for use by physicians trained and experienced in airway management.



Sizes and Corresponding Part Number

Part Number	Size	OD
V14750	14ch Vented 750mm Bougie	1.66mm

Technical details

Materials: High Density Polyethylene
Thermoplastic Elastomer (Adapter)

Latex status: Latex Free

Regulatory compliance

Meets the essential requirements of directive 93/42/EEC as amended 2007/47/EC for Class IIa Sterile medical Devices

Packaging: Unit pack - Single packed: film front fully peel-able pouch
Shelf pack - Cardboard tube of 10 pieces

Shelf Life: 5 Years from the date of sterilisation

Storage: There are no specific storage requirements

Manufacturing environment

Manufactured in a cleanroom constructed to comply with ISO 14644 Class 8 @ 0.5µ manned

Sterilisation

Terminally sterilised by Gamma Irradiation as per ISO 11137